

Description:

This Ultrasonic treatment helps to repair free radical damage, as well as reducing the loss of vital skin hydration making it essential for extremely dry, mature and photo-damaged skin. DermaSound Ultra works at the cellular level to provide maximum protection for the lipid membrane, averting cellular damage.

CLEANSE

1. Dispense 1 pump **Mega-Purifying Cleanser** to dry palms and apply to dry skin massaging for 1-2 minutes. Add water for slip, rinse with warm water.
2. Dispense 1 pump **Vitamin C Cleanser** to dry palms and apply to dry skin massaging for 1-2 minutes. Add water for slip, rinse with warm water.

TREAT

3. Mist skin with **Skin Recovery Mist**.
4. Dispense **OXYGEN Deep Pore Cleanser** and apply generously to the skin using a fan brush. Allow to remain on the skin for 3-5 minutes. Steam may be incorporated. Remove with warm towels.
5. Using the Derasound in Cavitation* Setting: KHz 90-99%, refer to DSU Cheat Sheet.
6. Dispense **Wrinkle Prescription Masque** and apply to the skin using a fan brush. Allow to remain on the skin for 10 minutes. Perform massage over masque will cold globes or quartz roller. Remove with cool towels.
7. Mist skin with **Skin Recovery Mist**.
8. Using the Derasound in Sonophoresis* setting, apply **Vitamin C Serum, Derma Pigment Skin Brightener, Skin Restoring Fulvic Elixir** and **Ultra Hydro Gel**.

BALANCE

9. Apply **Vitamin C Cream** and **Beauty Oil**. Dispense 1/4 tsp of each to palms and apply to face, neck, and décolleté.
10. Using the Derasound in Micro Amp* setting. Apply Ultra Hydro Gel if more slip is needed.
11. Apply a small, grain size amount of **CBD Regenerative Eye Cream** to entire eye area and massage until absorbed.

PROTECT

11. Finish with **Photo-Age Environmental Protection Gel 15**. Apply a small pea size amount to face, neck, décolleté and any other areas exposed to UVA/UVB Rays.

Recommended Treatment Time: 45 minutes